


Screen-Printed Carbon Electrodes

FMS-0008



Electrochemical Technique	Suitable
Chronoamperometry	<input checked="" type="checkbox"/>
Cyclic Voltammetry	<input checked="" type="checkbox"/>
Square Wave Voltammetry	<input checked="" type="checkbox"/>
Differential Pulse Voltammetry	<input checked="" type="checkbox"/>
Open Circuit Potentiometry	<input checked="" type="checkbox"/>
Electrochemical Impedance Spectroscopy	<input checked="" type="checkbox"/>
Others	Inquire

Sample Performance Data

Performance data

The following data describes typical expected performance from the sensors within a batch. Please contact the manufacturer for batch-to-batch performance information.

Method of analysis.

Analysis was performed by cyclic voltammetry using 5mM Potassium Ferri/ferrocyanide in 10mM PBS pH7.4

Settings

- Estart: 0.0V
- Evertex1: 1.0V
- Evertex2: -0.8V
- Steps: 0.01V
- Rate: 0.1V/s
- Cycles: 4

Results

The below voltammogram was obtained using the above settings. The initial scan (starting at 0.0V was removed)




Figure 1: 3 successive scans of mediator solution on sputtered carbon sensor

Sample Performance Data (cont.)

Between scan characteristics (within electrode)

Oxidation peak: 117.9uA at 0.41V,
Reduction peak: -101.7uA at -0.14V,

Peak repeatability (3 scans across 20 sensors within batch)

Oxidation peak coefficient of variation:

Stability 0.9% (scan to scan within sensor)
Repeatability 3.9% (between electrodes)

Reduction peak coefficient of variation:

Stability 0.2% (scan to scan within sensor)
Repeatability 3.4% (between electrodes)

Single use

The electrodes are designed to be single use, for best result use once then discard.

Technical Data

Parameter	Value
Working Electrode	Carbon
Counter Electrode	Carbon
Reference Electrode	Ag/AgCl
Substrate	PET
Conductive Tracks	Silver
Working Area	12.6 mm ²
Sample Volume	50 – 100 µL
Usage	Single use
Application	For R&D use only.

Engineering and usage data for FMS-0008 screen-printed carbon electrodes

For more information contact: info@FlexMedical-Solutions.com